

ARG66680 anti-Dengue virus Membrane (M) protein antibody [SQab19169]

Package: 100 µg
Store at: -20°C

Summary

Product Description	Mouse Monoclonal antibody [SQab19169] recognizes Dengue virus Membrane (M) protein
Tested Reactivity	DEN
Tested Application	ELISA, ICC/IF, WB
Specificity	This antibody recognizes DENV1, DENV2 and DENV3 Membrane protein
Host	Mouse
Clonality	Monoclonal
Clone	SQab19169
Isotype	IgG2b, kappa
Target Name	Dengue virus Membrane (M) protein
Species	Dengue virus
Immunogen	Dengue virus Membrane (M) protein.
Conjugation	Un-conjugated

Application Instructions

Application table	Application	Dilution
	ELISA	Assay-dependent
	ICC/IF	Assay-dependent
	WB	1 µg/ml
Application Note	ICC/IF: DENV3 not workable; ELISA: DENV1 and DENV3 not workable * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties

Form	Liquid
Purification	Purified by protein A/G.
Buffer	PBS (pH 7.4) and 0.01% Thimerosal.
Preservative	0.01% Thimerosal
Concentration	1 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.

Bioinformation

Highlight


Related products:

[Dengue virus antibodies](#); [Dengue virus ELISA Kits](#); [Dengue virus Duos / Panels](#);

Related news:


[Exploring Antiviral Immune Response](#)

Images



ARG66680 anti-Dengue virus Membrane (M) protein antibody [SQab19169] WB image

Western blot: 10 μ g of Vero and DENV2-infected Vero cell lysates stained with ARG66680 anti-Dengue virus Membrane (M) protein antibody [SQab19169] at 1 μ g/ml dilution.



ARG66680 anti-Dengue virus Membrane (M) protein antibody [SQab19169] WB image

Western blot: 10 μ g of DENV1-infected C636, DENV2-infected C636 and ZIKV-infected C636 cells (left to right). Cell lysates were stained with ARG66680 anti-Dengue virus Membrane (M) protein antibody [SQab19169] at 1 μ g/ml dilution.