

**ARG66777**  
**anti-RAB35 antibody**Package: 100 µl  
Store at: -20°C

### Summary

Product Description	Rabbit Polyclonal antibody recognizes RAB35
Tested Reactivity	Hu, Ms
Predict Reactivity	Bov, Rat
Tested Application	IHC-P, WB
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Target Name	RAB35
Species	Human
Immunogen	KLH-conjugated synthetic peptide within the center region of Human RAB35.
Conjugation	Un-conjugated
Alternate Names	Ras-related protein Rab-1C; RAB1C; Ras-related protein Rab-35; H-ray; GTP-binding protein RAY; RAY

### Application Instructions

Application table	Application	Dilution
	IHC-P	1:100
	WB	1:500 - 1:1000
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	
Positive Control	MCF7 and Raw264.7	
Observed Size	~ 23 kDa	


### Properties

Form	Liquid
Purification	Affinity purification with immunogen.
Buffer	0.42% Potassium phosphate (pH 7.3), 0.87% NaCl, 0.01% Sodium azide and 30% Glycerol.
Preservative	0.01% Sodium azide
Stabilizer	30% Glycerol
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.

Bioinformation


Gene Symbol	RAB35
Gene Full Name	RAB35, member RAS oncogene family
Function	<p>The small GTPases Rab are key regulators of intracellular membrane trafficking, from the formation of transport vesicles to their fusion with membranes. Rabs cycle between an inactive GDP-bound form and an active GTP-bound form that is able to recruit to membranes different sets of downstream effectors directly responsible for vesicle formation, movement, tethering and fusion. That Rab is involved in the process of endocytosis and is an essential rate-limiting regulator of the fast recycling pathway back to the plasma membrane. During cytokinesis, required for the postfurlrowing terminal steps, namely for intercellular bridge stability and abscission, possibly by controlling phosphatidylinositol 4,5-bis phosphate (PIP2) and SEPT2 localization at the intercellular bridge. May indirectly regulate neurite outgrowth. Together with TBC1D13 may be involved in regulation of insulin-induced glucose transporter SLC2A4/GLUT4 translocation to the plasma membrane in adipocytes. [UniProt]</p>
Calculated Mw	23 kDa
PTM	<p>AMPylation at Tyr-77 by L.pneumophila DrrA occurs in the switch 2 region and leads to moderate inactivation of the GTPase activity. It appears to prolong the lifetime of the GTP state of RAB1B by restricting access of GTPase effectors to switch 2 and blocking effector-stimulated GTP hydrolysis, thereby rendering RAB35 constitutively active.</p> <p>Phosphocholinated by L.pneumophila AnkX. Both GDP-bound and GTP-bound forms can be phosphocholinated. Phosphocholination inhibits the GEF activity of DENND1A. [UniProt]</p>
Cellular Localization	<p>Cell membrane; Lipid-anchor. Membrane, clathrin-coated pit. Cytoplasmic vesicle, clathrin-coated vesicle. Endosome. Melanosome. Note=Present on sorting endosomes and recycling endosome tubules. Tends to be enriched in PIP2-positive cell membrane domains. During mitosis, associated with the plasma membrane and present at the ingressing furrow during early cytokinesis as well as at the intercellular bridge later during cytokinesis. Identified in stage I to stage IV melanosomes. [UniProt]</p>

Images



ARG66777 anti-RAB35 antibody IHC-P image

Immunohistochemistry: Formalin-fixed and paraffin-embedded Human liver tissue stained with ARG66777 anti-RAB35 antibody.



#### ARG66777 anti-RAB35 antibody WB image

Western blot: MCF7 and Raw264.7 whole cell lysates stained with ARG66777 anti-RAB35 antibody.