

ARG70032 Human IL10 recombinant protein (Active) (His-tagged, C-ter)

Package: 100 µg, 20 µg
Store at: -20°C


Summary

Product Description	E. coli expressed, His-tagged (C-ter) Active Human IL10 recombinant protein
Tested Application	SDS-PAGE
Target Name	IL10
Species	Human
A.A. Sequence	Ser19 - Asn178
Expression System	E. coli
Activity	Active
Alternate Names	IL10A; GVHDS; TGIF; IL-10; Cytokine synthesis inhibitory factor; CSIF; Interleukin-10

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	PBS (pH 8.0)
Reconstitution	

Images



Human IL10 recombinant protein

ARG70032 Human IL10 recombinant protein (Active) (His-tagged, C-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70032 Human IL10 recombinant protein (Active) (His-tagged, C-ter).