

ARG70033**Human IL17 / IL17A recombinant protein (Active) (His-tagged, C-ter)**

Package: 100 µg, 20 µg

Store at: -20°C


Summary

Product Description	E. coli expressed, His-tagged (C-ter) Active Human IL17 / IL17A recombinant protein
Tested Application	SDS-PAGE
Target Name	IL17 / IL17A
Species	Human
A.A. Sequence	Ile20 - Ala155
Expression System	E. coli
Activity	Active
Activity Note	Determined by its ability to induce IL-6 secretion in 3T3 cells. The ED50 for this effect is < 6 ng/mL.
Alternate Names	IL-17A; Interleukin-17A; IL-17; CTLA8; IL17; CTLA-8; Cytotoxic T-lymphocyte-associated antigen 8

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	PBS (pH 8.0)
Reconstitution	

Images



ARG70033 Human IL17 / IL17A recombinant protein (Active) (His-tagged, C-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70033 Human IL17 / IL17A recombinant protein (Active) (His-tagged, C-ter).