

## ARG70041 Mouse TNF alpha recombinant protein (Active) (His-tagged, C-ter)

Package: 100 µg  
Store at: -20°C


### Summary

|                     |  |
|---------------------|--|
| Product Description | E. coli expressed, His-tagged (C-ter) Active Mouse TNF alpha recombinant protein   |
| Tested Application  | SDS-PAGE   |
| Target Name         | TNF alpha  |
| Species             | Mouse  |
| A.A. Sequence       | Leu80 - Leu235   |
| Expression System   | E. coli  |
| Activity            | Active   |
| Alternate Names     | Tumor necrosis factor ligand superfamily member 2; DIF; Cachectin; ICD2; ICD1; N-terminal fragment; TNF-a; TNFA; TNFSF2; TNF-alpha; Tumor necrosis factor; NTF |

### Properties

|                   |   |
|-------------------|---|
| Form              | Powder  |
| Purification Note | Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test. |
| Purity            | > 98% (by SDS-PAGE)   |
| Buffer            | PBS (pH 8.0)  |
| Reconstitution    |   |

### Images



ARG70041 Mouse TNF alpha recombinant protein (Active) (His-tagged, C-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70041 Mouse TNF alpha recombinant protein (Active) (His-tagged, C-ter).