

ARG70042 Human IL1 alpha recombinant protein (Active) (His-tagged, C-ter)

Package: 100 µg, 20 µg
Store at: -20°C


Summary

Product Description	E. coli expressed, His-tagged (C-ter) Active Human IL1 alpha recombinant protein
Tested Application	SDS-PAGE
Target Name	IL1 alpha
Species	Human
A.A. Sequence	Ser113 - Ala271
Expression System	E. coli
Activity	Active
Alternate Names	IL-1 alpha; Interleukin-1 alpha; IL1; IL1F1; Hematopoietin-1; IL1-ALPHA; IL-1A

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	PBS (pH 8.0)
Reconstitution	

Images



ARG70042 Human IL1 alpha recombinant protein (Active) (His-tagged, C-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70042 Human IL1 alpha recombinant protein (Active) (His-tagged, C-ter).