

## ARG70054 Human IL11 recombinant protein (Active) (His-tagged, N-ter)

Package: 100 µg, 20 µg  
Store at: -20°C


### Summary

Product Description	E. coli expressed, His-tagged (N-ter) Active Human IL11 recombinant protein
Tested Application	SDS-PAGE
Target Name	IL11
Species	Human
A.A. Sequence	Pro22 - Leu199
Expression System	E. coli
Activity	Active
Alternate Names	Oprelvekin; Interleukin-11; AGIF; Adipogenesis inhibitory factor; IL-11

### Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	PBS (pH 8.0)
Reconstitution	

### Images



ARG70054 Human IL11 recombinant protein (Active) (His-tagged, N-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70054 Human IL11 recombinant protein (Active) (His-tagged, N-ter).