

ARG70064 Human IL20 recombinant protein (Active) (His-tagged, C-ter)

Package: 100 µg, 20 µg
Store at: -20°C


Summary

Product Description	E. coli expressed, His-tagged (C-ter) Active Human IL20 recombinant protein
Tested Application	SDS-PAGE
Target Name	IL20
Species	Human
A.A. Sequence	Leu25 - Glu176
Expression System	E. coli
Activity	Active
Activity Note	Determined by its ability to chemoattract BaF3 cells transfected with human IL-20R alpha and IL-20R beta. The ED50 for this effect is < 0.2 ng/mL.
Alternate Names	Interleukin-20; IL-20; IL10D; Cytokine Zcyto10; ZCYTO10

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	PBS (pH 8.0)
Reconstitution	

Images



ARG70064 Human IL20 recombinant protein (Active) (His-tagged, C-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70064 Human IL20 recombinant protein (Active) (His-tagged, C-ter).