

ARG70069

Human IL25 / IL17E recombinant protein (Active) (His-tagged, C-ter)

Package: 100 µg, 20 µg

Store at: -20°C


Summary

Product Description	E. coli expressed, His-tagged (C-ter) Active Human IL25 / IL17E recombinant protein
Tested Application	SDS-PAGE
Target Name	IL25 / IL17E
Species	Human
A.A. Sequence	Tyr33 - Gly177
Expression System	E. coli
Activity	Active
Alternate Names	IL17E; Interleukin-17E; IL-25; Interleukin-25; IL-17E

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	PBS (pH 8.0)
Reconstitution	

Images



ARG70069 Human IL25 / IL17E recombinant protein (Active) (His-tagged, C-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70069 Human IL25 / IL17E recombinant protein (Active) (His-tagged, C-ter).