

**ARG70071****Human IL27B / EB13 recombinant protein (Active) (His-tagged, C-ter)**

Package: 100 µg, 20 µg

Store at: -20°C


### Summary

Product Description	E. coli expressed, His-tagged (C-ter) Active Human IL27B / EB13 recombinant protein
Tested Application	SDS-PAGE
Target Name	IL27B / EB13
Species	Human
A.A. Sequence	Arg21 - Lys229
Expression System	E. coli
Activity	Active
Activity Note	Determined by its ability to induce TF-1 cells proliferation using a concentration range of 20-200 ng/mL.
Alternate Names	EBV-induced gene 3 protein; IL-27 subunit beta; IL27B; Epstein-Barr virus-induced gene 3 protein; Interleukin-27 subunit beta; IL-27B

### Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 95% (by SDS-PAGE)
Buffer	PBS (pH 8.0)
Reconstitution	

### Images



Human IL27B / EB13 recombinant protein

ARG70071 Human IL27B / EB13 recombinant protein (Active) (His-tagged, C-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70071 Human IL27B / EB13 recombinant protein (Active) (His-tagged, C-ter).