

**ARG70072****Human IFN lambda 2 recombinant protein (Active) (His-tagged, C-ter)**

Package: 100 µg, 20 µg

Store at: -20°C


### Summary

Product Description	E. coli expressed, His-tagged (C-ter) Active Human IFN lambda 2 recombinant protein
Tested Application	SDS-PAGE
Target Name	IFN lambda 2
Species	Human
A.A. Sequence	Pro27 - Val200
Expression System	E. coli
Activity	Active
Activity Note	Determined by its ability to protect HepG2 cells infected with encephalomyocarditis (EMC) virus. The ED50 for this effect is < 20 ng/mL.
Alternate Names	Interleukin-28A; IFN-lambda-2; Interferon lambda-2; IL28A; IL-28A; Cytokine Zcyto20

### Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	PBS (pH 8.0)
Reconstitution	

### Images



Human IFN lambda 2 recombinant protein

ARG70072 Human IFN lambda 2 recombinant protein (Active) (His-tagged, C-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70072 Human IFN lambda 2 recombinant protein (Active) (His-tagged, C-ter).