

ARG70073 Human IL28B recombinant protein (Active) (His-tagged, C-ter)

Package: 100 µg, 20 µg
Store at: -20°C


Summary

Product Description	E. coli expressed, His-tagged (C-ter) Active Human IL28B recombinant protein
Tested Application	SDS-PAGE
Target Name	IL28B
Species	Human
A.A. Sequence	Val22 - Val196
Expression System	E. coli
Activity	Active
Activity Note	Determined by its ability to protect HepG2 cells infected with encephalomyocarditis (EMC) virus. The ED50 for this effect is < 5 ng/mL.
Alternate Names	Interleukin-28B; IFN-lambda-3; Interferon lambda-3; Cytokine Zcyto22; IL28B; IL28C; Interleukin-28C; IL-28B; IL-28C

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.01 EU/µg of the protein, as determined by the LAL test.
Purity	> 95% (by SDS-PAGE)
Buffer	PBS (pH 8.0)
Reconstitution	

Images



ARG70073 Human IL28B recombinant protein (Active) (His-tagged, C-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70073 Human IL28B recombinant protein (Active) (His-tagged, C-ter).