

**ARG70074****Human IL29 / IFN lambda recombinant protein (Active) (His-tagged, C-ter)**

Package: 100 µg, 20 µg

Store at: -20°C


## Summary

Product Description	E. coli expressed, His-tagged (C-ter) Active Human IL29 / IFN lambda recombinant protein
Tested Application	SDS-PAGE
Target Name	IL29 / IFN lambda
Species	Human
A.A. Sequence	Pro23 - Thr200
Expression System	E. coli
Activity	Active
Activity Note	Determined by its ability to induce IL-8 secretion in HuH7 cells. The ED50 for this effect is < 6 ng/mL.
Alternate Names	IFN-lambda-1; IL-29; Interferon lambda-1; Interleukin-29; Cytokine Zcyto21; IL29

## Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	PBS (pH 8.0)
Reconstitution	

## Images



Human IL29 / IFN lambda recombinant protein

ARG70074 Human IL29 / IFN lambda recombinant protein (Active) (His-tagged, C-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70074 Human IL29 / IFN lambda recombinant protein (Active) (His-tagged, C-ter).