

ARG70076 Human IL32 recombinant protein (Active) (His-tagged, C-ter)

Package: 100 µg, 20 µg
Store at: -20°C


Summary

Product Description	E. coli expressed, His-tagged (C-ter) Active Human IL32 recombinant protein
Tested Application	SDS-PAGE
Target Name	IL32
Species	Human
A.A. Sequence	Met1 - Lys131
Expression System	E. coli
Activity	Active
Activity Note	Determined by its ability to induce TNF alpha secretion in RAW264.7 cells. The ED50 for this effect is < 10 µg/mL.
Alternate Names	IL-32gamma; IL-32beta; IL-32delta; Interleukin-32; TAIa; TAIc; TAIb; NK4; TAI; IL-32; TAI; Natural killer cells protein 4; Tumor necrosis factor alpha-inducing factor; IL-32alpha

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	PBS (pH 8.0)
Reconstitution	

Images



Human IL32 (isoform alpha)
recombinant protein

ARG70076 Human IL32 recombinant protein (Active) (His-tagged, C-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70076 Human IL32 recombinant protein (Active) (His-tagged, C-ter).