

ARG70083

Human IL37 / IL1F7 / FIL1 zeta recombinant protein (Active) (His-tagged, C-ter)

Package: 100 µg, 20 µg


Store at: -20°C

Summary

Product Description	E. coli expressed, His-tagged (C-ter) Active Human IL37 / IL1F7 / FIL1 zeta recombinant protein
Tested Application	SDS-PAGE
Target Name	IL37 / IL1F7 / FIL1 zeta
Species	Human
A.A. Sequence	Lys53 - Asp218
Expression System	E. coli
Activity	Active
Activity Note	Determined by its ability to induce IL-8 secretion in human PBMCs. The ED50 for this effect is < 0.9 µg/mL.
Alternate Names	Interleukin-1 zeta; IL-1 zeta; FIL1; Interleukin-37; IL-1X; IL-1F7; FIL1 zeta; IL-1H4; Interleukin-1 homolog 4; FIL1Z; IL1H4; Interleukin-1 family member 7; Interleukin-23; IL-37; FIL1(ZETA); IL-1H; IL-1RP1; IL1RP1; IL1F7; Interleukin-1-related protein

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 95% (by SDS-PAGE)
Buffer	PBS (pH 8.0)
Reconstitution	



Human IL37 / IL1F7 / FIL1 zeta
recombinant protein

ARG70083 Human IL37 / IL1F7 / FIL1 zeta recombinant protein
(Active) (His-tagged, C-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70083 Human IL37 / IL1F7 / FIL1 zeta
recombinant protein (Active) (His-tagged, C-ter).