

ARG70086 Human BAFF recombinant protein (Active) (His-tagged, C-ter)

Package: 100 µg, 20 µg
Store at: -20°C


Summary

Product Description	E. coli expressed, His-tagged (C-ter) Active Human BAFF recombinant protein
Tested Application	SDS-PAGE
Target Name	BAFF
Species	Human
A.A. Sequence	Ala134 - Leu285
Expression System	E. coli
Activity	Active
Activity Note	Determined by its ability to induce IL-8 secretion in human PBMCs. The ED50 for this effect is < 0.5 ng/mL.
Alternate Names	BLYS; TALL1; Dendritic cell-derived TNF-like molecule; THANK; DTL; TNF- and APOL-related leukocyte expressed ligand 1; TNFSF20; CD257; BlyS; Tumor necrosis factor ligand superfamily member 13B; ZTNF4; CD antigen CD257; B lymphocyte stimulator; TALL-1; BAFF; B-cell-activating factor

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	PBS (pH 8.0)
Reconstitution	

Images



ARG70086 Human BAFF recombinant protein (Active) (His-tagged, C-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70086 Human BAFF recombinant protein (Active) (His-tagged, C-ter).