

## ARG70089 Human BMP4 recombinant protein (Active) (His-tagged, C-ter)

Package: 100 µg, 20 µg  
Store at: -20°C


### Summary

Product Description	E. coli expressed, His-tagged (C-ter) Active Human BMP4 recombinant protein
Tested Application	SDS-PAGE
Target Name	BMP4
Species	Human
A.A. Sequence	Lys303 - Arg408
Expression System	E. coli
Activity	Active
Activity Note	Determined by its ability to induce alkaline phosphatase production by ATDC5 cells. The ED50 for this effect is < 0.58 ng/mL.
Alternate Names	BMP2B; BMP-4; OFC11; Bone morphogenetic protein 2B; BMP2B1; Bone morphogenetic protein 4; MCOPS6; ZYME; BMP-2B

### Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 95% (by SDS-PAGE)
Buffer	20 mM sodium carbonate (pH 9.0)
Reconstitution	

### Images



ARG70089 Human BMP4 recombinant protein (Active) (His-tagged, C-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70089 Human BMP4 recombinant protein (Active) (His-tagged, C-ter).