

ARG70092 Human BMP7 recombinant protein (Active) (His-tagged, C-ter)

Package: 100 µg, 20 µg
Store at: -20°C


Summary

Product Description	E. coli expressed, His-tagged (C-ter) Active Human BMP7 recombinant protein
Tested Application	SDS-PAGE
Target Name	BMP7
Species	Human
A.A. Sequence	Met315 - His431
Expression System	E. coli
Activity	Active
Activity Note	Determined by its ability to induce alkaline phosphatase production by ATDC5 cells. The ED50 for this effect is < 0.65 µg/mL.
Alternate Names	BMP-7; OP-1; Osteogenic protein 1; Eptotermin alfa; Bone morphogenetic protein 7

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 95% (by SDS-PAGE)
Buffer	20 mM sodium citrate and 0.2 M NaCl (pH 3.5)
Reconstitution	

Images



ARG70092 Human BMP7 recombinant protein (Active) (His-tagged, C-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70092 Human BMP7 recombinant protein (Active) (His-tagged, C-ter).