

ARG70102 Human BMP16 recombinant protein (Active) (His-tagged, C-ter)

Package: 100 µg, 20 µg
Store at: -20°C


Summary

Product Description	E. coli expressed, His-tagged (C-ter) Active Human BMP16 recombinant protein
Tested Application	SDS-PAGE
Target Name	BMP16
Species	Human
A.A. Sequence	His238 - Leu347
Expression System	E. coli
Activity	Active
Activity Note	Determined by its ability to induce alkaline phosphatase production by ATDC5 cells. The ED50 for this effect is < 2.2 ng/mL.

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	20 mM sodium citrate and 0.2 M NaCl (pH 3.5)
Reconstitution	

Images



ARG70102 Human BMP16 recombinant protein (Active) (His-tagged, C-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70102 Human BMP16 recombinant protein (Active) (His-tagged, C-ter).