

ARG70104

Human FLT3 Ligand recombinant protein (Active) (His-tagged, C-ter)

Package: 100 µg, 20 µg

Store at: -20°C


Summary

| | |
|---------------------|--|
| Product Description | E. coli expressed, His-tagged (C-ter) Active Human FLT3 Ligand recombinant protein |
| Tested Application | SDS-PAGE |
| Target Name | FLT3 Ligand |
| Species | Human |
| A.A. Sequence | Thr27 - Ala181 |
| Expression System | E. coli |
| Activity | Active |
| Alternate Names | FLT3L; SL cytokine; Fms-related tyrosine kinase 3 ligand; Flt3L; Flt3 ligand; FL |

Properties

| | |
|-------------------|---|
| Form | Powder |
| Purification Note | Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test. |
| Purity | > 95% (by SDS-PAGE) |
| Buffer | PBS (pH 8.0) |
| Reconstitution | |

Images



Human FLT3 Ligand recombinant protein

ARG70104 Human FLT3 Ligand recombinant protein (Active) (His-tagged, C-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70104 Human FLT3 Ligand recombinant protein (Active) (His-tagged, C-ter).