

ARG70107 Human TGF beta 2 recombinant protein (Active) (His-tagged, C-ter)

Package: 100 µg, 20 µg
Store at: -20°C


Summary

Product Description	E. coli expressed, His-tagged (C-ter) Active Human TGF beta 2 recombinant protein
Tested Application	SDS-PAGE
Target Name	TGF beta 2
Species	Human
A.A. Sequence	Ala303 - Ser414
Expression System	E. coli
Activity	Active
Alternate Names	TGF-beta2; Cetermin; LDS4; Polyergin; Glioblastoma-derived T-cell suppressor factor; G-TSF; Transforming growth factor beta-2; LAP; TGF-beta-2; BSC-1 cell growth inhibitor

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.01 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	20 mM Sodium citrate and 0.2 M NaCl (pH 3.5).
Reconstitution	

Images



ARG70107 Human TGF beta 2 recombinant protein (Active) (His-tagged, C-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70107 Human TGF beta 2 recombinant protein (Active) (His-tagged, C-ter).