

ARG70108 Human FGF acidic recombinant protein (Active) (His-tagged, C-ter)

Package: 100 µg, 20 µg
Store at: -20°C


Summary

Product Description	E. coli expressed, His-tagged (C-ter) Active Human FGF acidic recombinant protein
Tested Application	SDS-PAGE
Target Name	FGF acidic
Species	Human
A.A. Sequence	Phe16 - Asp155
Expression System	E. coli
Activity	Active
Alternate Names	HBGF1; FGF-1; FGF-alpha; GLIO703; ECGFA; ECGFB; ECGF; Endothelial cell growth factor; FGFA; Acidic fibroblast growth factor; Heparin-binding growth factor 1; Fibroblast growth factor 1; HBGF-1; AFGF; aFGF; ECGF-beta

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	PBS (pH 8.0)
Reconstitution	

Images



ARG70108 Human FGF acidic recombinant protein (Active) (His-tagged, C-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70108 Human FGF acidic recombinant protein (Active) (His-tagged, C-ter).