

ARG70109 Human FGF basic recombinant protein (Active) (His-tagged, N-ter)

Package: 100 µg, 20 µg
Store at: -20°C

Summary

Product Description	E. coli expressed, His-tagged (N-ter) Active Human FGF basic recombinant protein
Tested Application	SDS-PAGE
Target Name	FGF basic
Species	Human
A.A. Sequence	Ala135 - Ser288
Expression System	E. coli
Activity	Active
Alternate Names	FGF-2; Fibroblast growth factor 2; bFGF; FGFB; Heparin-binding growth factor 2; BFGF; HBGF-2; Basic fibroblast growth factor

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	PBS containing 0.01% sarkosyl (pH 8.0)
Reconstitution	

Images



Human FGF basic recombinant protein

ARG70109 Human FGF basic recombinant protein (Active) (His-tagged, N-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70109 Human FGF basic recombinant protein (Active) (His-tagged, N-ter).