

ARG70111 Human FGF4 recombinant protein (Active) (His-tagged, C-ter)

Package: 100 µg, 20 µg
Store at: -20°C


Summary

Product Description	E. coli expressed, His-tagged (C-ter) Active Human FGF4 recombinant protein
Tested Application	SDS-PAGE
Target Name	FGF4
Species	Human
A.A. Sequence	Gly25 - Leu206
Expression System	E. coli
Activity	Active
Alternate Names	K-FGF; HST-1; FGF-4; Transforming protein KS3; HSTF-1; Heparin secretory-transforming protein 1; Fibroblast growth factor 4; KFGF; HST; Heparin-binding growth factor 4; HBGF-4; HSTF1

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 95% (by SDS-PAGE)
Buffer	PBS containing 0.1% sarkosyl (pH 8.0)
Reconstitution	

Images



ARG70111 Human FGF4 recombinant protein (Active) (His-tagged, C-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70111 Human FGF4 recombinant protein (Active) (His-tagged, C-ter).