

## ARG70112 Human FGF5 recombinant protein (Active) (His-tagged, C-ter)

Package: 100 µg, 20 µg  
Store at: -20°C


### Summary

Product Description	E. coli expressed, His-tagged (C-ter) Active Human FGF5 recombinant protein
Tested Application	SDS-PAGE
Target Name	FGF5
Species	Human
A.A. Sequence	Ala18 - Gly268
Expression System	E. coli
Activity	Active
Alternate Names	FGF-5; TCMGLY; Fibroblast growth factor 5; Heparin-binding growth factor 5; Smag-82; HBGF-5

### Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 95% (by SDS-PAGE)
Buffer	PBS (pH 8.0)
Reconstitution	

### Images



ARG70112 Human FGF5 recombinant protein (Active) (His-tagged, C-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70112 Human FGF5 recombinant protein (Active) (His-tagged, C-ter).