

ARG70114

Human FGF8 (isoform A) recombinant protein (Active) (His-tagged, C-ter)

Package: 100 µg, 20 µg

Store at: -20°C


Summary

Product Description	E. coli expressed, His-tagged (C-ter) Active Human FGF8 (isoform A) recombinant protein
Tested Application	SDS-PAGE
Target Name	FGF8 (isoform A)
Species	Human
A.A. Sequence	Gln52 - Arg223
Expression System	E. coli
Activity	Active
Activity Note	Determined by its ability to induce 3T3 cells proliferation. The ED50 for this effect is < 2 µg/mL.
Alternate Names	Fibroblast growth factor 8; HH6; AIGF; FGF-8; HBGF-8; KAL6; Heparin-binding growth factor 8; Androgen-induced growth factor

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	PBS (pH 8.0)
Reconstitution	

Images



Human FGF8 (isoform A)
recombinant protein

ARG70114 Human FGF8 (isoform A) recombinant protein (Active)
(His-tagged, C-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70114 Human FGF8 (isoform A)
recombinant protein (Active) (His-tagged, C-ter).