

ARG70115 Human FGF9 recombinant protein (Active) (His-tagged, C-ter)

Package: 100 µg, 20 µg
Store at: -20°C


Summary

Product Description	E. coli expressed, His-tagged (C-ter) Active Human FGF9 recombinant protein
Tested Application	SDS-PAGE
Target Name	FGF9
Species	Human
A.A. Sequence	Pro3 - Ser208
Expression System	E. coli
Activity	Active
Activity Note	Determined by its ability to induce 3T3 cells proliferation. The ED50 for this effect is < 2 ng/mL.
Alternate Names	Fibroblast growth factor 9; Glia-activating factor; FGF-9; HBFG-9; HBGF-9; GAF; Heparin-binding growth factor 9; SYNS3

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 95% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Reconstitution	

Images



ARG70115 Human FGF9 recombinant protein (Active) (His-tagged, C-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70115 Human FGF9 recombinant protein (Active) (His-tagged, C-ter).