

ARG70117

Human FGF11 (isoform 1) recombinant protein (Active) (His-tagged, C-ter)

Package: 100 µg, 20 µg

Store at: -20°C


Summary

Product Description	E. coli expressed, His-tagged (C-ter) Active Human FGF11 (isoform 1) recombinant protein
Tested Application	SDS-PAGE
Target Name	FGF11 (isoform 1)
Species	Human
A.A. Sequence	Met1 - Pro225
Expression System	E. coli
Activity	Active
Activity Note	Determined by its ability to induce 3T3 cells proliferation. The ED50 for this effect is < 0.2 ng/mL.
Alternate Names	FHF-3; Fibroblast growth factor homologous factor 3; Fibroblast growth factor 11; FGF-11; FHF3

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.01 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Reconstitution	

Images



Human FGF11 (isoform 1)
recombinant protein

ARG70117 Human FGF11 (isoform 1) recombinant protein (Active)
(His-tagged, C-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70117 Human FGF11 (isoform 1)
recombinant protein (Active) (His-tagged, C-ter).