

## ARG70120 Human FGF13 recombinant protein (Active) (His-tagged, C-ter)

Package: 100 µg, 20 µg  
Store at: -20°C


### Summary

Product Description	E. coli expressed, His-tagged (C-ter) Active Human FGF13 recombinant protein
Tested Application	SDS-PAGE
Target Name	FGF13
Species	Human
A.A. Sequence	Met1 - Thr245
Expression System	E. coli
Activity	Active
Activity Note	Determined by its ability to induce 3T3 cells proliferation. The ED50 for this effect is < 160 ng/mL.
Alternate Names	FHF-2; FHF2; Fibroblast growth factor 13; Fibroblast growth factor homologous factor 2; FGF-13; FGF2

### Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Reconstitution	

### Images



ARG70120 Human FGF13 recombinant protein (Active) (His-tagged, C-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70120 Human FGF13 recombinant protein (Active) (His-tagged, C-ter).