

ARG70126 Human FGF20 recombinant protein (Active) (His-tagged, C-ter)

Package: 100 µg, 20 µg
Store at: -20°C


Summary

| | |
|---------------------|---|
| Product Description | E. coli expressed, His-tagged (C-ter) Active Human FGF20 recombinant protein |
| Tested Application | SDS-PAGE |
| Target Name | FGF20 |
| Species | Human |
| A.A. Sequence | Pro3 - Thr211 |
| Expression System | E. coli |
| Activity | Active |
| Activity Note | Determined by its ability to induce 3T3 cells proliferation. The ED50 for this effect is 1.3-3.2 ng/mL. The specific activity of recombinant human FGF-20 is $> 2 \times 10^5$ IU/mg. |
| Alternate Names | Fibroblast growth factor 20; FGF-20; RHDA2 |

Properties

| | |
|-------------------|---|
| Form | Powder |
| Purification Note | Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test. |
| Purity | $> 98\%$ (by SDS-PAGE) |
| Buffer | PBS (pH 8.0) |
| Reconstitution | |

Images



ARG70126 Human FGF20 recombinant protein (Active) (His-tagged, C-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70126 Human FGF20 recombinant protein (Active) (His-tagged, C-ter).