

ARG70128**Human FGF22 recombinant protein (Active) (His-SUMO tagged, N-ter)**

Package: 100 µg, 20 µg

Store at: -20°C


Summary

Product Description	E. coli expressed, His-SUMO tagged (N-ter) Active Human FGF22 recombinant protein
Tested Application	SDS-PAGE
Target Name	FGF22
Species	Human
A.A. Sequence	Thr23 - Ser170
Expression System	E. coli
Activity	Active
Activity Note	Determined by its ability to induce 3T3 cells proliferation. The ED50 for this effect is < 2 ng/mL.
Alternate Names	Fibroblast growth factor 22; FGF-22

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	PBS (pH 8.0)
Reconstitution	

Images



ARG70128 Human FGF22 recombinant protein (Active) (His-SUMO tagged, N-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70128 Human FGF22 recombinant protein (Active) (His-SUMO tagged, N-ter).