

ARG70139

Human Tpo / Thrombopoietin recombinant protein (His-tagged, N-ter)

Package: 100 µg, 20 µg

Store at: -20°C

Summary

Product Description	E. coli expressed, His-tagged (N-ter) Human Tpo / Thrombopoietin recombinant protein
Tested Application	SDS-PAGE
Target Name	Tpo / Thrombopoietin
Species	Human
A.A. Sequence	Ser22 - Leu195
Expression System	E. coli
Alternate Names	C-mpl ligand; Thrombopoietin; MPLLG; TPO; ML; Megakaryocyte colony-stimulating factor; Megakaryocyte growth and development factor; MGDF; Myeloproliferative leukemia virus oncogene ligand; MKCSF; THCYT1

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.01 EU/µg of the protein, as determined by the LAL test.
Purity	> 95% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Reconstitution	

Images



Human Tpo / Thrombopoietin recombinant protein

ARG70139 Human Tpo / Thrombopoietin recombinant protein (His-tagged, N-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70139 Human Tpo / Thrombopoietin recombinant protein (His-tagged, N-ter).