

ARG70144 Human Neurturin recombinant protein (Active) (His-tagged, C-ter)

Package: 100 µg, 20 µg
Store at: -20°C


Summary

| | |
|---------------------|---|
| Product Description | E. coli expressed, His-tagged (C-ter) Active Human Neurturin recombinant protein |
| Tested Application | SDS-PAGE |
| Target Name | Neurturin |
| Species | Human |
| A.A. Sequence | Ala96 - Val197 |
| Expression System | E. coli |
| Activity | Active |
| Activity Note | Determined by its ability to induce proliferation in SH-SY5Y cells. The ED50 for this effect is < 50 ng/mL. |
| Alternate Names | Neurturin; NTN |

Properties

| | |
|-------------------|---|
| Form | Powder |
| Purification Note | Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test. |
| Purity | > 98% (by SDS-PAGE) |
| Buffer | PBS (pH 7.4) |
| Reconstitution | |

Images



ARG70144 Human Neurturin recombinant protein (Active) (His-tagged, C-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70144 Human Neurturin recombinant protein (Active) (His-tagged, C-ter).