

ARG70147 Human Galectin 3 recombinant protein (Active) (His-tagged, N-ter)

Package: 100 µg, 20 µg
Store at: -20°C


Summary

Product Description	E. coli expressed, His-tagged (N-ter) Active Human Galectin 3 recombinant protein
Tested Application	SDS-PAGE
Target Name	Galectin 3
Species	Human
A.A. Sequence	Ala2 - Ile250
Expression System	E. coli
Activity	Active
Activity Note	Determined by its ability to chemoattract human PBMC using a concentration range of 2.5-25 µg/mL. Note: Results may vary from different PBMC donors.
Alternate Names	Laminin-binding protein; Gal-3; L-31; GALBP; Galactoside-binding protein; MAC2; GAL3; GALIG; Mac-2 antigen; CBP 35; Galectin-3; CBP35; Galactose-specific lectin 3; IgE-binding protein; L31; 35 kDa lectin; Carbohydrate-binding protein 35; Lectin L-29

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 95% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Reconstitution	

Images



Human Galectin 3
recombinant protein

ARG70147 Human Galectin 3 recombinant protein (Active) (His-tagged, N-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70147 Human Galectin 3 recombinant protein (Active) (His-tagged, N-ter).