

ARG70148 Human Galectin 4 recombinant protein (Active) (His-tagged, N-ter)

Package: 100 µg, 20 µg
Store at: -20°C


Summary

Product Description	E. coli expressed, His-tagged (N-ter) Active Human Galectin 4 recombinant protein
Tested Application	SDS-PAGE
Target Name	Galectin 4
Species	Human
A.A. Sequence	Ala2 - Ile323
Expression System	E. coli
Activity	Active
Activity Note	Determined by its ability to agglutinate human red blood cells. The ED50 for this effect is < 8 µg/mL.
Alternate Names	GAL4; L36LBP; L-36 lactose-binding protein; Gal-4; Galectin-4; Lactose-binding lectin 4; Antigen NY-CO-27

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Reconstitution	

Images



Human Galectin 4
recombinant protein

ARG70148 Human Galectin 4 recombinant protein (Active) (His-tagged, N-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70148 Human Galectin 4 recombinant protein (Active) (His-tagged, N-ter).