

## ARG70149 Human Galectin 7 recombinant protein (Active) (His-tagged, N-ter)

Package: 100 µg, 20 µg  
Store at: -20°C


### Summary

Product Description	E. coli expressed, His-tagged (N-ter) Active Human Galectin 7 recombinant protein
Tested Application	SDS-PAGE
Target Name	Galectin 7
Species	Human
A.A. Sequence	Ser2 - Phe136
Expression System	E. coli
Activity	Active
Activity Note	Determined by its ability to agglutinate human red blood cells. The ED50 for this effect is < 2 µg/mL
Alternate Names	GAL7; Gal-7; HKL-14; LGALS7A; Galectin-7; p53-induced gene 1 protein; PI7

### Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Reconstitution	

### Images



ARG70149 Human Galectin 7 recombinant protein (Active) (His-tagged, N-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70149 Human Galectin 7 recombinant protein (Active) (His-tagged, N-ter).