

ARG70153 Mouse IL1RA recombinant protein (Active) (His-tagged, C-ter)

Package: 100 µg, 20 µg
Store at: -20°C


Summary

Product Description	E. coli expressed, His-tagged (C-ter) Active Mouse IL1RA recombinant protein
Tested Application	SDS-PAGE
Target Name	IL1RA
Species	Mouse
A.A. Sequence	Arg27 - Gln178
Expression System	E. coli
Activity	Active
Activity Note	Determined by its ability to inhibit IL-1 alpha -dependent proliferation in D10.G4.1 cells. The ED50 for this effect is < 50 ng/mL.
Alternate Names	IL1 inhibitor; Interleukin-1 receptor antagonist protein; ICIL-1RA; MVCD4; IL1RA; DIRA; IL-1RN; Anakinra; IL1F3; IRAP; IL-1ra; IL-1ra3

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Reconstitution	

Images



ARG70153 Mouse IL1RA recombinant protein (Active) (His-tagged, C-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70153 Mouse IL1RA recombinant protein (Active) (His-tagged, C-ter).