

## ARG70155 Mouse IL3 recombinant protein (Active) (His-tagged, C-ter)

Package: 100 µg, 20 µg  
Store at: -20°C


### Summary

Product Description	E. coli expressed, His-tagged (C-ter) Active Mouse IL3 recombinant protein
Tested Application	SDS-PAGE
Target Name	IL3
Species	Mouse
A.A. Sequence	Asp33 - Cys166
Expression System	E. coli
Activity	Active
Alternate Names	MCGF; Mast cell growth factor; P-cell-stimulating factor; Hematopoietic growth factor; IL-3; Interleukin-3; Multipotential colony-stimulating factor; MULTI-CSF

### Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 95% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Reconstitution	

### Images



ARG70155 Mouse IL3 recombinant protein (Active) (His-tagged, C-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70155 Mouse IL3 recombinant protein (Active) (His-tagged, C-ter).