

ARG70160 Mouse IL9 recombinant protein (Active) (His-tagged, N-ter)

Package: 100 µg, 20 µg
Store at: -20°C


Summary

Product Description	E. coli expressed, His-tagged (N-ter) Active Mouse IL9 recombinant protein
Tested Application	SDS-PAGE
Target Name	IL9
Species	Mouse
A.A. Sequence	Gln19 - Pro144
Expression System	E. coli
Activity	Active
Alternate Names	IL-9; HP40; Interleukin-9; T-cell growth factor P40; Cytokine P40; P40

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Reconstitution	

Images



ARG70160 Mouse IL9 recombinant protein (Active) (His-tagged, N-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70160 Mouse IL9 recombinant protein (Active) (His-tagged, N-ter).