

ARG70163
Mouse IL12A / IL12 p35 recombinant protein (Active) (His-tagged, C-ter)

Package: 100 µg, 20 µg


Store at: -20°C

Summary

Product Description	E. coli expressed, His-tagged (C-ter) Active Mouse IL12A / IL12 p35 recombinant protein
Tested Application	SDS-PAGE
Target Name	IL12A / IL12 p35
Species	Mouse
A.A. Sequence	Arg23 - Ala215
Expression System	E. coli
Activity	Active
Activity Note	Determined by its ability to induce proliferation in T-cell enriched PBMC. The ED50 for this effect is < 0.2 ng/mL.
Alternate Names	CLMF; Cytotoxic lymphocyte maturation factor 35 kDa subunit; NFSK; NK cell stimulatory factor chain 1; CLMF p35; NKSF1; IL-12 subunit p35; IL-12A; Interleukin-12 subunit alpha; P35

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Reconstitution	



ARG70163 Mouse IL12A / IL12 p35 recombinant protein (Active) (His-tagged, C-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70163 Mouse IL12A / IL12 p35 recombinant protein (Active) (His-tagged, C-ter).