

## ARG70164 Mouse IL12B / IL12 p40 recombinant protein (Active) (His-tagged, C-ter)


Package: 100 µg, 20 µg  
Store at: -20°C

### Summary

Product Description	E. coli expressed, His-tagged (C-ter) Active Mouse IL12B / IL12 p40 recombinant protein
Tested Application	SDS-PAGE
Target Name	IL12B / IL12 p40
Species	Mouse
A.A. Sequence	Met23 - Ser335
Expression System	E. coli
Activity	Active
Activity Note	Determined by its ability to induce proliferation in T-cell enriched PBMC. The ED50 for this effect is < 0.3 ng/mL.
Alternate Names	CLMF; Interleukin-12 subunit beta; NK cell stimulatory factor chain 2; NKSF; CLMF p40; NKSF2; Cytotoxic lymphocyte maturation factor 40 kDa subunit; CLMF2; IL-12 subunit p40; IMD29; IMD28; IL-12B

### Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Reconstitution	



Mouse IL12B / IL12 p40  
recombinant protein

ARG70164 Mouse IL12B / IL12 p40 recombinant protein (Active) (His-tagged, C-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70164 Mouse IL12B / IL12 p40 recombinant protein (Active) (His-tagged, C-ter).