

ARG70167 Mouse IL17 / IL17A recombinant protein (Active) (His-tagged, C-ter)

Package: 100 µg, 20 µg
Store at: -20°C


Summary

Product Description	E. coli expressed, His-tagged (C-ter) Active Mouse IL17 / IL17A recombinant protein
Tested Application	SDS-PAGE
Target Name	IL17 / IL17A
Species	Mouse
A.A. Sequence	Ala26 - Ala158
Expression System	E. coli
Activity	Active
Alternate Names	IL-17A; Interleukin-17A; IL-17; CTLA8; IL17; CTLA-8; Cytotoxic T-lymphocyte-associated antigen 8

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	20 mM sodium citrate and 0.2 M NaCl (pH 4.5)
Reconstitution	

Images



ARG70167 Mouse IL17 / IL17A recombinant protein (Active) (His-tagged, C-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70167 Mouse IL17 / IL17A recombinant protein (Active) (His-tagged, C-ter).