

ARG70168 Mouse IL17F recombinant protein (Active) (His-tagged, C-ter)

Package: 100 µg, 20 µg
Store at: -20°C


Summary

Product Description	E. coli expressed, His-tagged (C-ter) Active Mouse IL17F recombinant protein
Tested Application	SDS-PAGE
Target Name	IL17F
Species	Mouse
A.A. Sequence	Arg29 - Ala161
Expression System	E. coli
Activity	Active
Activity Note	Determined by its ability to induce IL-6 secretion in 3T3 cells. The ED50 for this effect is < 100 ng/mL.
Alternate Names	Interleukin-17F; IL-17F; CANDF6; ML-1; ML1; Cytokine ML-1

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	20 mM sodium citrate (pH 4.5)
Reconstitution	

Images



ARG70168 Mouse IL17F recombinant protein (Active) (His-tagged, C-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70168 Mouse IL17F recombinant protein (Active) (His-tagged, C-ter).