

## ARG70170 Mouse IL19 recombinant protein (Active) (His-tagged, C-ter)

Package: 100 µg, 20 µg  
Store at: -20°C


### Summary

Product Description	E. coli expressed, His-tagged (C-ter) Active Mouse IL19 recombinant protein
Tested Application	SDS-PAGE
Target Name	IL19
Species	Mouse
A.A. Sequence	Leu25 - Ala176
Expression System	E. coli
Activity	Active
Activity Note	Determined by its ability to induce proliferation in BaF3 cells transfected with human IL-20 R alpha and human IL-20 R beta. The ED50 for this effect is < 0.6 ng/mL.
Alternate Names	Interleukin-19; IL-10C; ZMDA1; IL-19; MDA1; NG.1; Melanoma differentiation-associated protein-like protein

### Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Reconstitution	

### Images



ARG70170 Mouse IL19 recombinant protein (Active) (His-tagged, C-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70170 Mouse IL19 recombinant protein (Active) (His-tagged, C-ter).