

## ARG70175 Mouse IL25 / IL17E recombinant protein (Active) (His-tagged, C-ter)

Package: 100 µg, 20 µg  
Store at: -20°C


### Summary

|                     |                                                                                                          |
|---------------------|----------------------------------------------------------------------------------------------------------|
| Product Description | E. coli expressed, His-tagged (C-ter) Active Mouse IL25 / IL17E recombinant protein                      |
| Tested Application  | SDS-PAGE                                                                                                 |
| Target Name         | IL25 / IL17E                                                                                             |
| Species             | Mouse                                                                                                    |
| A.A. Sequence       | Val17 - Ala169                                                                                           |
| Expression System   | E. coli                                                                                                  |
| Activity            | Active                                                                                                   |
| Activity Note       | Determined by its ability to induce CXCL1 secretion in HT29 cells. The ED50 for this effect is < 1 ng/mL |
| Alternate Names     | IL17E; Interleukin-17E; IL-25; Interleukin-25; IL-17E                                                    |

### Properties

|                   |                                                                               |
|-------------------|-------------------------------------------------------------------------------|
| Form              | Powder                                                                        |
| Purification Note | Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test. |
| Purity            | > 98% (by SDS-PAGE)                                                           |
| Buffer            | 20 mM sodium citrate and 0.2 M NaCl (pH 4.5)                                  |
| Reconstitution    |                                                                               |

### Images



ARG70175 Mouse IL25 / IL17E recombinant protein (Active) (His-tagged, C-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70175 Mouse IL25 / IL17E recombinant protein (Active) (His-tagged, C-ter).