

ARG70179 Mouse IL27 recombinant protein (Active) (His-tagged, C-ter)

Package: 100 µg, 20 µg
Store at: -20°C


Summary

Product Description	E. coli expressed, His-tagged (C-ter) Active Mouse IL27 recombinant protein
Tested Application	SDS-PAGE
Target Name	IL27
Species	Mouse
A.A. Sequence	Phe29 - Ser234
Expression System	E. coli
Activity	Active
Activity Note	Determined by its ability to protect HepG2 cells infected with encephalomyocarditis (EMC) virus. The ED50 for this effect is < 5 ng/mL.
Alternate Names	IL27; Interleukin 27; IL27A; P28; IL27p28; IL-27A; IL-27; IL30; Interleukin-27 Subunit Alpha; IL-27 Subunit Alpha; Interleukin-30; MGC71873; IL-27-A; IL27-A; IL-27 P28 Subunit; Interleukin 30

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Reconstitution	

Images



ARG70179 Mouse IL27 recombinant protein (Active) (His-tagged, C-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70179 Mouse IL27 recombinant protein (Active) (His-tagged, C-ter).