

ARG70180 Mouse IL33 recombinant protein (Active) (His-tagged, C-ter)

Package: 100 µg, 20 µg
Store at: -20°C


Summary

Product Description	E. coli expressed, His-tagged (C-ter) Active Mouse IL33 recombinant protein
Tested Application	SDS-PAGE
Target Name	IL33
Species	Mouse
A.A. Sequence	Ser109 - Ile266
Expression System	E. coli
Activity	Active
Alternate Names	95-270; NF-HEV; Interleukin-33; C9orf26; IL1F11; 99-270; Interleukin-1 family member 11; IL-33; IL-1F11; Nuclear factor from high endothelial venules; NFEHEV; DVS27; 109-270

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Reconstitution	

Images



ARG70180 Mouse IL33 recombinant protein (Active) (His-tagged, C-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70180 Mouse IL33 recombinant protein (Active) (His-tagged, C-ter).