

ARG70187

Mouse FLT3 Ligand recombinant protein (Active) (His-tagged, C-ter)

Package: 100 µg, 20 µg

Store at: -20°C


Summary

Product Description	E. coli expressed, His-tagged (C-ter) Active Mouse FLT3 Ligand recombinant protein
Tested Application	SDS-PAGE
Target Name	FLT3 Ligand
Species	Mouse
A.A. Sequence	Gly27 - Arg188
Expression System	E. coli
Activity	Active
Activity Note	Determined by its ability to induce proliferation in BaF3 cells transfected with mouse Flt-3. The ED50 for this effect is < 2 ng/mL.
Alternate Names	FLT3L; SL cytokine; Fms-related tyrosine kinase 3 ligand; Flt3L; Flt3 ligand; FL

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Reconstitution	

Images



ARG70187 Mouse FLT3 Ligand recombinant protein (Active) (His-tagged, C-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70187 Mouse FLT3 Ligand recombinant protein (Active) (His-tagged, C-ter).