

## ARG70189 Mouse FGF basic recombinant protein (Active) (His-tagged, N-ter)

Package: 100 µg, 20 µg  
Store at: -20°C


### Summary

Product Description	E. coli expressed, His-tagged (N-ter) Active Mouse FGF basic recombinant protein
Tested Application	SDS-PAGE
Target Name	FGF basic
Species	Mouse
A.A. Sequence	Ala11 - Ser154
Expression System	E. coli
Activity	Active
Alternate Names	FGF-2; Fibroblast growth factor 2; bFGF; FGFB; Heparin-binding growth factor 2; BFGF; HBGF-2; Basic fibroblast growth factor

### Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	PBS containing 0.01% sarkosyl (pH 8.0)
Reconstitution	

### Images



ARG70189 Mouse FGF basic recombinant protein (Active) (His-tagged, N-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70189 Mouse FGF basic recombinant protein (Active) (His-tagged, N-ter).