

ARG70190 Mouse IGF1 recombinant protein (Active) (His-tagged, C-ter)

Package: 100 µg, 20 µg
Store at: -20°C

Summary

Product Description	E. coli expressed, His-tagged (C-ter) Active Mouse IGF1 recombinant protein
Tested Application	SDS-PAGE
Target Name	IGF1
Species	Mouse
A.A. Sequence	Gly49 - Ala118
Expression System	E. coli
Activity	Active
Alternate Names	MGF; Insulin-like growth factor I; Mechano growth factor; Somatomedin-C; IGF1; IGF-I

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	PBS (pH 8.0)
Reconstitution	

Images



Mouse IGF1
recombinant protein

ARG70190 Mouse IGF1 recombinant protein (Active) (His-tagged, C-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70190 Mouse IGF1 recombinant protein (Active) (His-tagged, C-ter).